




Redipor[®] BetaBag

Redipor[®] BetaBag, from AnalytiChem UK (formerly Cherwell Laboratories), supports continuous manufacturing by reducing time and risk during product transfer.





redipor

an analytichem brand

This product works in conjunction with the Getinge La Calhène DPTE® Alpha Port of 190 mm diameter as part of an isolator or RABS.

The Redipor® BetaBag is an addition to the well-established Redipor® range of environmental monitoring solutions, ensuring that the highest levels of production and patient safety is achieved.

This technology is designed as a solution to support continuous manufacturing cycles, focusing on:

- ▶ Contamination-free transfer of Redipor® plated media to maintain rapid production
- ▶ Fully validated and guaranteed transfer process
- ▶ Significantly reduced risk of biological and particulate contamination
- ▶ Save on validated and on-site cleaning times





Optimized productivity

Irradiated Contact or Settle Plates

The transfer system enables you to gain access to a game-changing solution, that optimizes cleanroom operations to ensure the integrity of sterile products, which is critical to your environment.

Our development enables you to stay ahead in an increasingly competitive market, demonstrating your commitment to:

- ▶ Quality
- ▶ Safety
- ▶ Efficiency, which is in optimizing productivity in pharmaceutical production

Improve your Environmental Monitoring

The purpose of the Redipor® BetaBag is to ensure the safe transfer of ready-to-use gamma-irradiated prepared media, into your sterile manufacturing environment, while reducing risk, cost, and time.

The advantages of the Redipor® BetaBag also include:

- ▶ Can be attached to an alpha port up to four times without compromising its integrity
- ▶ Convenient storage
- ▶ Additional space within an isolator or RABS
- ▶ Eliminating requirements for vaporised hydrogen peroxide (VHP) decontamination

This method truly simplifies your workflow giving the process greater adaptability, while ensuring maximum safety of your product, and patients.

Comparison of the decontamination process using AnalytiChem UK's VHP compatible barrier pack product against the new Redipor® BetaBag when obtaining plates for environmental monitoring (EM).

Current Four Stage Process

Introduction of plates into isolator

- ▶ Plates are removed from their outer packaging, and put into the isolator via the transfer hatch; then another layer of packaging is removed.

Decontamination cycle

- ▶ The decontamination cycle is performed according to the customer's SOP.

Storage, usage and disposal

- ▶ The plates are removed from their inner packaging and are stored within the isolator, where they're used for environmental monitoring. Any unused plates are disposed of.

Repeat process

- ▶ If any additional plates are needed, the decontamination process must be performed again.

Three Stage Process using a Redipor® BetaBag

(Eliminating Decontamination Stage)

Introduction of plates into isolator

- ▶ The BetaBag is attached to the alpha port on the isolator

Storage

- ▶ Plates are stored in the bag

Continuous EM

- ▶ Each sleeve of plates is removed from the bag as needed, enabling continuous environmental monitoring




What is in the bag?

AnalytiChem UK provides a comprehensive range of plates specifically designed for environmental monitoring in isolators and critical cleanrooms.

- ▶ Contact or Settle Plates for maximum safety, reliability and convenience
- ▶ Prepared according to European and US pharmacopoeia recommendations
- ▶ Available with a GS1-compliant data matrix barcode
- ▶ Produced in cleanrooms and gamma irradiated in the final packaging
- ▶ Long shelf life
- ▶ Available with 4 neutralizers for inactivation of a wide range of disinfectants and inhibiting VHP residue
- ▶ Room temperature stability





Ordering Information

The following products are available from AnalytiChem UK (formerly Cherwell Laboratories).

If you require a different formulation, fill volume or packaging format, contact **sales.uk@analytichem.com** to discuss your unique requirements.

For availability outside the UK, visit our website or contact **info@analytichem.com** to learn which Redipor BetaBag products are offered in your country.


Part Number	Description	Qty / Pack
5.0818	Irradiated Tryptone Soya Agar 90mm Petri dish single vent clear 18ml fill. Flow wrap + label. Sleeve of 10	1 Redipor® BetaBag (10 x 10)
5.0842	Irradiated Tryptone Soya 1.6% Agar + Neutralizer No.4 90mm Petri dish single vent clear data matrix 18ml fill. Barrier pack. Sleeve of 10	1 Redipor® BetaBag (10 x 10)
5.0819	Irradiated Tryptone Soya 1.6% Agar + Neutralizer No.4 55mm Contact plate triple vent clear 17ml fill. Flow wrap + label. Sleeve of 10	1 Redipor® BetaBag (10 x 10)
5.0809	Irradiated Tryptone Soya 2% Agar + Neutralizer No.4 90mm Petri dish single vent clear expiry date 27ml fill. Flow wrap + label. Sleeve of 10	1 Redipor® BetaBag (10 x 10)
5.0810	Irradiated Tryptone Soya 2% Agar + Neutralizer No.4 55mm Contact plate triple vent clear expiry date 17ml fill. Flow wrap + label. Sleeve of 10	1 Redipor® BetaBag (10 x 10)



Terms & Conditions

We reserve the right to amend specifications without notice. E&OE





redipor 
an analyti**chem** brand

AnalytiChem UK Ltd.
sales.uk@analytichem.com
+44 (0) 1869 355500

7 & 8 Launton Business Centre,
Murdock Road, Bicester,
OX26 4XB, UK



© 2025 AnalytiChem. All rights reserved.
All trademarks and registered trademarks
are the property of their respective owners.
www.analytichem.com