

Universal Operator Broth Transfer Kit

For USP<797>, ISO 13408-1, PIC/S, GMP Guidelines

Designed based on the UK NHS-protocol, the Universal Operator Broth Transfer kit is ideal for the verification of the aseptic handling techniques of your personnel at your sterile compounding facility or other clean space.

This kit includes:

- ▶ 1 X TSB glass break ampoule
- ▶ 1 X TSB, 50 mL vial
- ▶ 1 X TBS, 100 mL infusion bag
- ▶ 3 X Sterile vials

All culture media in the kit undergoes robust quality testing according to EP/USP, including analysis of growth promotion, pH and sterility. The product is accompanied by a certificate of analysis so you can confidently assess your operator's ability to maintain the sterility of materials during handling.


Description	Part Number
Universal Operator Broth Transfer Kit 1	UOBTk1




UNIVERSAL OPERATOR BROTH TRANSFER TEST

Schematic Representation of Universal Broth Test Procedure:


Transfer from bag to syringes




Transfer from vial to bag



Transfer from ampoule to vial



Transfer from bag to vial



Transfer from vial to vial

