



QbigSeT

Sterile Packaging Solutions

QbigSeT is a versatile customizable packaging solution developed for sterile media preparation. Each is treated with Ethylene Oxide (EO) to ensure the sterility of the external surfaces of bottles or tubes sealed with screw or crimp caps.

Sizes & Specifications

Box Size	Dimensions (L × W × H)	Bottle Formats	Max Batch Size (Sets)
Size 1	22.0 × 18.0 × 10.5 cm	20 × Alpha bottle 90 mL, screw cap	60
		25 × Alpha bottle 60 mL, screw cap	48
Size 2	25.0 × 20.0 × 11.5 cm	20 × Alpha bottle 125 mL, screw cap	60
		20 × Infusion bottle 100 mL, fels cap	60
		120 × Tube 15 mL, screw cap	22
Size 3	32.5 × 26.0 × 15.5 cm	20 × Alpha bottle 300 mL, screw cap	35 to 40
Size 5	32.5 × 24.5 × 18.0 cm	12 × Alpha bottle 500 mL, screw cap	30 to 35
		12 × Infusion bottle 500 mL, fels cap	30 to 35

Customization Options

Can be tailored to meet customer-specific requirements:

- ▶ Bottle or tube type
- ▶ Media formulation and volume
- ▶ Cap type
- ▶ Number of sterile barriers

Product Codes

Size	Code	Description
Size 1	KQST1E1	Box for 20 bottles of 90 mL, EO treated, 1 bag
	KQST1E2	Box for 20 bottles of 90 mL, EO treated, 2 bags
Size 2	KQST2E1	Box for 20 bottles of 100 or 125 mL or 120 tubes of 15 mL, EO treated, 1 bag
	KQST2E2	Box for 20 bottles of 100 or 125 mL or 120 tubes of 15 mL, EO treated, 2 bags
Size 3	KQST3E1	Box for 20 bottles of 300 mL, EO treated, 1 bag
	KQST3E2	Box for 20 bottles, EO treated, 2 bags
Size 5	KQST5E1	Box for 12 bottles, EO treated, 1 bag
	KQST5E2	Box for 12 bottles, EO treated, 2 bags

Ordering Examples

Article code	Article description	Packaging
KQST1E2-11	Size 1 number 11 (10 x 45 mL TSB + 1% Tween and 10 x 45 mL FTM + 1% Tween)	Size 1 , double wrapped, with: 10 x K235F045FZ and 10 x K236F045FY
KQST2E1-11	Size 2 number 11 (20 x 100 mL Rinsing fluid D)	Size 2 , single wrapped, with: 20 x K177F100GM
KQST2E1-13	Size 2 number 13 (10 x 100 mL TSB and 10 x 100 mL FTM for filtration)	Size 2 , single wrapped, with: 10 x K111F100CU and 10 x K899F100CT
KQST2E1-15	Size 2 number 15 (20 x 100 mL TSB)	Size 2 , single wrapped, with: 20 x K111F100CU
KQST2E1-17	Size 2 number 17 (20 x 100 mL TSB)	Size 2 , single wrapped, with: 20 x K111F100GM
KQST2E1-18	Size 2 number 18 (20 x 100 mL FTM for filtration)	Size 2 , single wrapped, with: 20 x K899F100GL
KQST2E1-19	Size 2 number 19 (20 x 100 mL FTM for filtration)	Size 2 , single wrapped, with: 20 x K899F100CT
KQST2E1-22	Size 2 number 22 (20 x 100 mL Rinsing fluid D)	Size 2 , single wrapped, with: 20 x K177F100CU
KQST2E2-14	Size 2 number 14 (20 x 100 mL TSB)	Size 2 , double wrapped, with: 20 x K111F100GM
KQST2E2-16	Size 2 number 16 (20 x 100 mL FTM)	Size 2 , double wrapped, with: 20 x K119F100GL
KQST2E2-19	Size 2 number 19 (10 x 100 mL TSB and 10 x 100 mL FTM)	Size 2 , double wrapped, with: 10 x K111F100GM and 10 x K119F100GL
KQST2E2-20	Size 2 number 20 (20 x 100 mL Dest. water)	Size 2 , double wrapped, with: 20 x K180F100CU
KQST2E2-27	Size 2 number 27 (10 x 100 mL TSB and 10 x 100 mL FTM for filtration)	Size 2 , double wrapped, with: 10 x K111F100CU and 10 x K899F100CT
KQST3E2-14	Size 3 number 14 (20 x 100 mL TSB)	Size 3 , double wrapped, with: 20 x K111F100GZ
KQST5E1-14	Size 5 number 14 (12 x 300 mL GPP)	Size 5 , single wrapped, with: 12 x K287F300HQ
KQST5E2-12	Size 5 number 12 (12 x 400 mL Rinsing fluid A)	Size 5 , double wrapped, with: 12 x K100F400DR
KQST5E2-14	Size 5 number 14 (12 x 300 mL Rinsing fluid A)	Size 5 , double wrapped, with: 12 x K100F300DR
KQST5E2-15	Size 5 number 15 (12 x 100 mL Rinsing fluid D)	Size 5 , double wrapped, with: 12 x K177F100DR
KQST5E2-16	Size 5 number 16 (12 x 300 mL Rinsing fluid D)	Size 5 , double wrapped, with: 12 x K177F300DR
KQST5E2-22	Size 5 number 22 (12 x 300 mL Rinsing fluid K)	Size 5 , double wrapped, with: 12 x K348F300DR
KQST5E3-11	Size 5 number 11 (12 x 100 mL Rinsing fluid A)	Size 5 , triple wrapped, with: 12 x K100F300DR
KQST5E3-12	Size 5 number 12 (12 x 300 mL Rinsing fluid K)	Size 5 , triple wrapped, with: 12 x K348F300DR

