


Peristaltic Pump Tubing



Contents

2-Stop Options:

PVC	3
Santoprene	4
SolventFlex	5
Viton	6
Silicon	6

3-Stop Options

PVC	7
Santoprene	8
SolventFlex	9
Viton	.10
Silicon	.10

Our portfolio includes a broad range of peristaltic pump tubing to cover most ICP-OES and ICP-MS applications. Our pump tubing is quality controlled with strict dimensional tolerances, ensuring consistent flow rate from tube to tube. Explore a selection of our most popular tubes below.

PVC 2-Stop Options 15 tubes / pack		Length (mm) / Interstop distance (mm)					
		288mm/72mm		455mm/150 mm		900mm/151 mm	320mm/95 mm
Stop Colors	I.D.	Not Flared	Flared	Not Flared	Flared	Not Flared	Flared
Orange/Red	0.19 mm	023-030-002	023-130-002	022-030-002	022-130-002		025-130-002
Orange/Blue	0.25 mm	023-030-003	023-130-003	022-030-003	022-130-003	022-030-603	025-130-003
Orange/Green	0.38 mm	023-030-004	023-130-004	022-030-004	022-130-004	022-030-604	025-130-004
Green/Yellow	0.44 mm	023-030-024	023-130-024	022-030-024	022-130-024		025-130-124
Orange/Yellow	0.51 mm	023-030-005	023-130-005	022-030-005	022-130-005	022-030-605	025-130-005
Orange/White	0.64 mm	023-030-006	023-130-006	022-030-006	022-130-006	022-030-606	025-130-006
Black/Black	0.76 mm	023-030-007	023-130-007	022-030-007	022-130-007	022-030-607	025-130-007
Orange/Orange	0.89 mm	023-030-008	023-130-008	022-030-008	022-130-008	022-030-608	
White/Black	0.95 mm	023-030-023	023-130-023	022-030-023	022-130-023	022-030-623	
White/White	1.02 mm	023-030-009	023-130-009	022-030-009	022-130-009	022-030-609	
Red/Red	1.14 mm	023-030-010	023-130-010	022-030-010	022-130-010	022-030-610	
Gray/Gray	1.30 mm	023-030-011	023-130-011	022-030-011	022-130-011	022-030-611	
Yellow/Yellow	1.42 mm	023-030-012	023-130-012	022-030-012	022-130-012	022-030-612	
Yellow/Blue	1.52 mm	023-030-019	023-130-019	022-030-019	022-130-019	022-030-619	
Blue/Blue	1.65 mm	023-030-013	023-130-013	022-030-013	022-130-013	022-030-613	
Blue/Green	1.75 mm	023-030-022	023-130-022	022-030-022	022-130-022		
Green/Green	1.85 mm	023-030-014	023-130-014	022-030-014	022-130-014	022-030-614	
Purple/Purple	2.06 mm	023-030-015	023-130-015	022-030-015	022-130-015	022-030-615	
Purple/Black	2.29 mm	023-030-016	023-130-016	022-030-016			
Purple/Orange	2.54 mm	023-030-017	023-130-017	022-030-017	022-130-017		
Purple/White	2.79 mm	023-030-018	023-130-018	022-030-018		022-030-618	
Black/White	3.18 mm	023-030-021	023-130-021	022-030-021	022-130-021	022-030-621	

* Length is 900 mm.

Santoprene 2-Stop Options 15 tubes / pack		Length (mm) / Interstop distance (mm)	
		288mm/72 mm	400mm/150 mm
Stop Colors	I.D.	Not Flared	
Orange/Red	0.19 mm	023-030-402	022-030-402
Orange/Blue	0.25 mm	023-030-403	022-030-403
Orange/Green	0.38 mm	023-030-404	022-030-404
Green/Yellow	0.44 mm	023-030-424	022-030-424
Orange/Yellow	0.51 mm	023-030-405	022-030-405
Orange/White	0.64 mm	023-030-406	022-030-406
Black/Black	0.76 mm	023-030-407	022-030-407
Orange/Orange	0.89 mm	023-030-408	022-030-408
White/Black	0.95 mm	023-030-423	022-030-423
White/White	1.02 mm	023-030-409	022-030-409
Red/Red	1.14 mm	023-030-410	022-030-410
Gray/Gray	1.30 mm	023-030-411	022-030-411
Yellow/Yellow	1.42 mm	023-030-412	022-030-412
Yellow/Blue	1.52 mm	023-030-419	022-030-419
Blue/Blue	1.65 mm	023-030-413	022-030-413
Blue/Green	1.75 mm	023-030-422	022-030-422
Green/Green	1.85 mm	023-030-414	022-030-414
Purple/Purple	2.06 mm	023-030-415	022-030-415
Purple/Black	2.29 mm	023-030-416	022-030-416
Purple/Orange	2.54 mm	023-030-417	022-030-417
Purple/White	2.79 mm	023-030-418	022-030-418
Black/White	3.18 mm	023-030-421	022-030-421

SolventFlex 2-Stop Options 15 tubes / pack		Length (mm) / Interstop distance (mm)					
		288mm/72 mm		455mm/140 mm		310mm/165mm	320mm/95 mm
Stop Colors	I.D.	Not Flared	Flared	Not Flared	Flared	Not Flared	Flared
Orange/Red	0.19 mm	023-030-102	023-130-102	022-030-102	022-130-102		025-130-102
Orange/Blue	0.25 mm	023-030-103	023-130-103	022-030-103	022-130-103		
Orange/Green	0.38 mm	023-030-104	023-130-104	022-030-104	022-130-104		025-130-104
Green/Yellow	0.44 mm	023-030-124	023-130-124	022-030-124	022-130-124		025-130-124
Orange/Yellow	0.51 mm	023-030-105	023-130-105	022-030-105	022-130-105	022-037-005	
Orange/White	0.64 mm	023-030-106	023-130-106	022-030-106	022-130-106		025-130-106
Black/Black	0.76 mm	023-030-107	023-130-107	022-030-107	022-130-107	022-037-007	025-130-107
Orange/Orange	0.89 mm	023-030-108	023-130-108	022-030-108	022-130-108		
White/Black	0.95 mm	023-030-123	023-130-123	022-030-123	022-130-123		
White/White	1.02 mm	023-030-109	023-130-109	022-030-109	022-130-109		
Red/Red	1.14 mm	023-030-110	023-130-110	022-030-110	022-130-110	022-037-010	025-130-110
Gray/Gray	1.30 mm	023-030-111	023-130-111	022-030-111	022-130-111		
Yellow/Yellow	1.42 mm	023-030-112	023-130-112	022-030-112	022-130-112		
Yellow/Blue	1.52 mm	023-030-119	023-130-119	022-030-119	022-130-119		
Blue/Blue	1.65 mm	023-030-113	023-130-113	022-030-113	022-130-113		
Blue/Green	1.75 mm	023-030-122	023-130-122	022-030-122	022-130-122		
Green/Green	1.85 mm	023-030-114	023-130-114	022-030-114	022-130-114		
Purple/Purple	2.06 mm	023-030-115	023-130-115	022-030-115	022-130-115		
Purple/Black	2.29 mm	023-030-116	023-130-116	022-030-116	022-130-116		
Purple/Orange	2.54 mm	023-030-117	023-130-117	022-030-117	022-130-117		
Purple/White	2.79 mm	023-030-118	023-130-118	022-030-118	022-130-118		022-030-118
Black/White	3.18 mm	023-030-121	023-130-121	022-030-121	022-130-121		

Viton 2-Stop Options 12 tubes / pack		Length (mm) / Interstop distance (mm)	
		175mm/127 mm	455mm/140 mm
Stop Colors	I.D.	Not Flared	
Orange/Yellow	0.51 mm	020-030-305	020-031-305
Orange/White	0.64 mm	020-030-306	020-031-306
Black/Black	0.76 mm	020-030-307	020-031-307
Orange/Orange	0.89 mm	020-030-308	020-031-308
White/White	1.02 mm	020-030-309	020-031-309
Red/Red	1.14 mm	020-030-310	020-031-310
Gray/Gray	1.30 mm	020-030-311	020-031-311
Yellow/Yellow	1.42 mm	020-030-312	020-031-312
Yellow/Blue	1.52 mm	020-030-319	
Blue/Blue	1.65 mm	020-030-313	
Green/Green	1.85 mm	020-030-314	
Purple/Purple	2.06 mm	020-030-315	
Purple/Black	2.29 mm	020-030-316	
Purple/White	2.79 mm	020-030-318	
Yellow/Yellow	1.40 mm	020-030-352	

Silicon 2-Stop Options 6 tubes / pack		Length (mm) / Interstop distance (mm)	
		455mm/140 mm	
Stop Colors	I.D.	Not Flared	
Orange/Green	0.38 mm	020-030-204	
Orange/Yellow	0.51 mm	020-030-205	
Orange/White	0.64 mm	020-030-206	
Black/Black	0.76 mm	020-030-207	
Orange/Orange	0.89 mm	020-030-208	
White/White	1.02 mm	020-030-209	
Red/Red	1.14 mm	020-030-210	
Gray/Gray	1.30 mm	020-030-211	
Yellow/Yellow	1.42 mm	020-030-212	
Yellow/Blue	1.52 mm	020-030-219	
Blue/Blue	1.65 mm	020-030-213	
Green/Green	1.85 mm	020-030-214	
Purple/Purple	2.06 mm	020-030-215	
Purple/Black	2.29 mm	020-030-216	
Purple/White	2.79 mm	020-030-218	

PVC 3-Stop Options 15 tubes/pack		Length (mm) / Interstop distance (mm)			
		395 mm/95 mm		455 mm/72 mm	
Stop Colors	I.D.	Not Flared	Flared	Not Flared	Flared
Orange/Red/Orange	0.19mm	024-033-002	024-133-002	022-033-002	022-133-002
Orange/Blue/Orange	0.25mm	024-033-003	024-133-003	022-033-003	022-133-003
Orange/Green/Orange	0.38mm	024-033-004	024-133-004	022-033-004	022-133-004
Green/Yellow/Green	0.44mm	024-033-024	024-133-024		
Orange/Yellow/Orange	0.51mm	024-033-005	024-133-005	022-033-005	022-133-005
Orange/White/Orange	0.64mm	024-033-006	024-133-006	022-033-006	022-133-006
Black/Black/Black	0.76mm	024-033-007	024-133-007	022-033-007	022-133-007
Orange/Orange/Orange	0.89mm	024-033-008	024-133-008	022-033-008	022-133-008
White/Black/White	0.95mm	024-033-023	024-133-023	022-033-023	022-133-023
White/White/White	1.02mm	024-033-009	024-133-009	022-033-009	022-133-009
Red/Red/Red	1.14mm	024-033-010	024-133-010	022-033-010	022-133-010
Gray/Gray/Gray	1.30mm	024-033-011	024-133-011	022-033-011	022-133-011
Yellow/Yellow/Yellow	1.42mm	024-033-012	024-133-012	022-033-012	
Yellow/Blue/Yellow	1.52mm	024-033-019	024-133-019	022-033-019	022-133-019
Blue/Blue/Blue	1.65mm	024-033-013	024-133-013	022-033-013	
Blue/Green/Blue	1.75mm	024-033-022	024-133-022		
Green/Green/Green	1.85mm	024-033-014	024-133-014	022-033-014	
Purple/Purple/Purple	2.06mm	024-033-015	024-133-015	022-033-015	
Purple/Black/Purple	2.29mm	024-033-016	024-133-016	022-033-016	
Purple/Orange/Purple	2.54mm	024-033-017	024-133-017	022-033-017	022-133-017
Purple/White/Purple	2.79mm	024-033-018	024-133-018	022-033-018	022-133-018
Black/White/Black	3.18mm	024-033-021	024-133-021	022-033-021	

Santoprene 3-Stop Options 15 tubes/pack		Length (mm) / Interstop distance (mm)	
		395 mm/95 mm	400 mm/72 mm
Stop Colors	I.D.	Not Flared	
Orange/Red/Orange	0.19mm	024-033-402	022-033-402
Orange/Blue/Orange	0.25mm	024-033-403	022-033-403
Orange/Green/Orange	0.38mm	024-033-404	022-033-404
Green/Yellow/Green	0.44mm	024-033-424	022-033-424
Orange/Yellow/Orange	0.51mm	024-033-405	022-033-405
Orange/White/Orange	0.64mm	024-033-406	022-033-406
Black/Black/Black	0.76mm	024-033-407	022-033-407
Orange/Orange/Orange	0.89mm	024-033-408	022-033-408
White/Black/White	0.95mm	024-033-423	022-033-423
White/White/White	1.02mm	024-033-409	022-033-409
Red/Red/Red	1.14mm	024-033-410	022-033-410
Gray/Gray/Gray	1.30mm	024-033-411	022-033-411
Yellow/Yellow/Yellow	1.42mm	024-033-412	022-033-412
Yellow/Blue/Yellow	1.52mm	024-033-419	022-033-419
Blue/Blue/Blue	1.65mm	024-033-413	022-033-413
Blue/Green/Blue	1.75mm	024-033-422	022-033-422
Green/Green/Green	1.85mm	024-033-414	022-033-414
Purple/Purple/Purple	2.06mm	024-033-415	022-033-415
Purple/Black/Purple	2.29mm	024-033-416	022-033-416
Purple/Orange/Purple	2.54mm	024-033-417	022-033-417
Purple/White/Purple	2.79mm	024-033-418	022-033-418
Black/White/Black	3.18mm	024-033-421	022-033-421

SolventFlex 3-Stop Options 15 tubes/pack		Length (mm) / Interstop distance (mm)				
		395mm/95 mm		455mm/82 mm		
Stop Colors	I.D.	Not Flared	Flared	Not Flared	Flared	
Orange/Red/Orange	0.19mm	024-033-102	024-133-102	022-033-102	022-133-102	
Orange/Blue/Orange	0.25mm	024-033-103	024-133-103	022-033-103	022-133-103	
Orange/Green/Orange	0.38mm	024-033-104	024-133-104	022-033-104	022-133-104	
Green/Yellow/Green	0.44mm	024-033-124	024-133-124	022-033-124	022-133-124	
Orange/Yellow/Orange	0.51mm	024-033-105	024-133-105	022-033-105	022-133-105	
Orange/White/Orange	0.64mm	024-033-106	024-133-106	022-033-106	022-133-106	
Black/Black/Black	0.76mm	024-033-107	024-133-107	022-033-107	022-133-107	
Orange/Orange/Orange	0.89mm	024-033-108	024-133-108	022-033-108	022-133-108	
White/Black/White	0.95mm	024-033-123	024-133-123	022-033-123	022-133-123	
White/White/White	1.02mm	024-033-109	024-133-109	022-033-109	022-133-109	
Red/Red/Red	1.14mm	024-033-110	024-133-110	022-033-110	022-133-110	
Gray/Gray/Gray	1.30mm	024-033-111	024-133-111	022-033-111	022-133-111	
Yellow/Yellow/Yellow	1.42mm	024-033-112	024-133-112	022-033-112	022-133-112	
Yellow/Blue/Yellow	1.52mm	024-033-119	024-133-119	022-033-119	022-133-119	
Blue/Blue/Blue	1.65mm	024-033-113	024-133-113	022-033-113	022-133-113	
Blue/Green/Blue	1.75mm	024-033-122	024-133-122	022-033-122	022-133-122	
Green/Green/Green	1.85mm	024-033-114	024-133-114	022-033-114	022-133-114	
Purple/Purple/Purple	2.06mm	024-033-115	024-133-115	022-033-115	022-133-115	
Purple/Black/Purple	2.29mm	024-033-116	024-133-116	022-033-116	022-133-116	
Purple/Orange/Purple	2.54mm	024-033-117	024-133-117	022-033-117	022-133-117	
Purple/White/Purple	2.79mm	024-033-118	024-133-118	022-033-118	022-133-118	
Black/White/Black	3.18mm	024-033-121	024-133-121	022-033-121	022-133-121	

Viton 3-Stop Options 12 tubes/pack		Length (mm) / Interstop distance (mm)
		400mm/82 mm
Stop Colors	I.D.	Not Flared
Orange/Yellow/Orange	0.51 mm	020-033-305
Orange/White/Orange	0.64mm	020-033-306
Black/Black/Black	0.76 mm	020-033-307
Orange/Orange/Orange	0.89 mm	020-033-308
White/White/White	1.02 mm	020-033-309
Red/Red/Red	1.14 mm	020-033-310
Yellow/Yellow/Yellow	1.42 mm	020-033-312
Yellow/Blue/Yellow	1.52 mm	020-033-319
Blue/Blue/Blue	1.65 mm	020-033-313
Green/Green/Green	1.85 mm	020-033-314
Purple/Purple/Purple	2.06 mm	020-033-315
Purple/Black/Purple	2.29mm	020-033-316
Purple/Orange/Purple	2.54mm	020-033-317
Purple/White/Purple	2.79 mm	020-033-318

Silicon 3-Stop Options 6 tubes/pack		Length (mm) / Interstop distance (mm)
		400mm/67 mm
Stop Colors	I.D.	Not Flared
Orange/White/Orange	0.64 mm	020-033-206
Black/Black/Black	0.76 mm	020-033-207
Orange/Orange/Orange	0.89 mm	020-033-208
White/White/White	1.02 mm	020-033-209
Red/Red/Red	1.14 mm	020-033-210
Yellow/Yellow/Yellow	1.42 mm	020-033-212
Yellow/Blue/Yellow	1.52 mm	020-033-219
Blue/Blue/Blue	1.65 mm	020-033-213
Green/Green/Green	1.85 mm	020-033-214
Purple/Purple/Purple	2.06 mm	020-033-215
Purple/Black/Purple	2.29mm	020-033-216
Purple/White/Purple	2.79 mm	020-033-218



Americas & Asia
salesna@analytichem.com

EMEA
info@analytichem.com

Australia
sales@ore.com.au

Additional Info & Contacts
www.analytichem.com

© 2025 AnalytiChem. All rights reserved. All trademarks and registered trademarks are the property of their respective owners.