

DigiTUBEs

Disposable digestion tubes

DigiTUBEs are disposable digestion tubes suitable for use with commercially available graphite digestion block systems. They are manufactured with ICP-MS Grade Polypropylene and specially designed with a flat bottom to provide more direct contact with the base of digestion block wells, resulting in a faster heat-up time and more efficient sample temperature response times.

Every lot of **DigiTUBE** digestion tubes receives a Certificate of Volume analysis confirming the product's compliance with Class A Specification according to ASTM E1272 and ISO 384 per our ISO 17025 accredited procedure.



Additionally, **DigiTUBE**s are now available with a Trace Metals Certificate, providing verified results for over 65 elements to guarantee ultra-low contamination and the compliance required for critical trace metals analyses.

For clean applications, **DigiTUBE**s are also available in shrink-wrapped bundles to protect against dust infiltration during storage and handling.

Key Features

- ▶ Ultra-low leachable metal content to eliminate contamination risk
- ▶ Optional RackLock feature enables one-handed closure of screw caps
- ▶ Class A graduations for easy and accurate dilution
- ▶ Disposable, eliminating the need to clean the vessel
- ▶ Leak-free screw cap available in variety of colors



*DigiTUBE*s can replace 5 vessels in the lab - digestion beaker, graduated cylinder, volumetric flask, autosampler tube, and storage container. Save time and lower costs with the use of *DigiTUBE*s.

Bulk *DigiTUBE*s

Description	Quantity	Part Number
15 mL with RackLock, blue caps	1620	010-515-261
15 mL with RackLock, no caps	1620	010-515-262
15 mL without RackLock, blue caps	1620	010-515-263
15 mL without RackLock, no caps	1620	010-515-264
15 mL with RackLock, blue caps	540	010-515-361
15 mL without RackLock, blue caps	540	010-515-363
50 mL with RackLock, blue caps	750	010-500-261
50 mL with RackLock, clear caps	750	010-500-265
50 mL with RackLock, no caps	750	010-500-262
50 mL without RackLock, blue caps	750	010-500-263
50 mL without RackLock, clear caps	750	010-500-266
50 mL without RackLock, no caps	750	010-500-264
100 mL with RackLock, blue caps	300	010-501-261
100 mL with RackLock, no caps	300	010-501-262
100 mL without RackLock, blue caps	300	010-501-263
100 mL without RackLock, no caps	300	010-501-264

Bulk Caps

Screw Caps for	Description	Quantity	Part Number
15 mL <i>DigiTUBE</i> s	Screw Cap, Blue for 15 mL <i>DigiTUBE</i> s	540	010-515-060
	Screw Cap, Clear for 15 mL <i>DigiTUBE</i> s	540	010-515-172
50 mL <i>DigiTUBE</i> s	Blue Caps for 50 mL <i>DigiTUBE</i> s	250	010-500-060
	Yellow Caps for 50 mL <i>DigiTUBE</i> s	250	010-500-130
	Orange Caps for 50 mL <i>DigiTUBE</i> s	250	010-500-140
	Red Caps for 50 mL <i>DigiTUBE</i> s	250	010-500-150
	Clear Caps for 50 mL <i>DigiTUBE</i> s	250	010-500-172
100 mL <i>DigiTUBE</i> s	Blue Caps for 100 mL <i>DigiTUBE</i> s	100	010-501-060
	Clear Caps for 100 mL <i>DigiTUBE</i> s	100	010-501-172

HandyPAKs

HandyPAKs include multiple racks of *DigiTUBE*s (50 mL or 100 mL), shrink wrapped to eliminate contamination during storage or transport. Caps are provided in a separate bag. The **HandyPAK Plus** option contains pre-capped *DigiTUBE*s.

	Description	Quantity	Part Number
HandyPAK	50 mL <i>DigiTUBE</i> s with RackLock, blue caps	30 x 25 tubes	010-500-362
	50 mL <i>DigiTUBE</i> s without RackLock, blue caps	30 x 25 tubes	010-500-364
	50 mL <i>DigiTUBE</i> s with RackLock, clear caps	30 x 25 tubes	010-500-366
	100 mL <i>DigiTUBE</i> s with RackLock, blue caps	9 x 34 tubes	010-501-362
	100 mL <i>DigiTUBE</i> s without RackLock, blue caps	9 X 34 tubes	010-501-364
HandyPAK Plus	50 mL <i>DigiTUBE</i> s with RackLock, blue caps	20 x 24 tubes	010-500-461
	50 mL <i>DigiTUBE</i> s without RackLock, blue caps	20 x 24 tubes	010-500-463
	100 mL <i>DigiTUBE</i> s with RackLock, blue caps	24 x 12 tubes	010-501-462
	100 mL <i>DigiTUBE</i> s without RackLock, blue caps	24 x 12 tubes	010-501-464
HandySleeve	50mL <i>DigiTUBE</i> s without RackLock, blue cap	30 x 25 tubes	010-500-563
	50mL <i>DigiTUBE</i> s without RackLock, clear cap	30 x 25 tubes	010-500-564
	50mL <i>DigiTUBE</i> s without RackLock, red cap	30 x 25 tubes	010-500-565
	50mL <i>DigiTUBE</i> s with RackLock, blue cap	30 x 25 tubes	010-500-561
	50mL <i>DigiTUBE</i> s with RackLock, clear cap	30 x 25 tubes	010-500-562
	100mL <i>DigiTUBE</i> s without RackLock, blue cap	30x10 tubes	010-501-563
	100mL <i>DigiTUBE</i> s without RackLock, clear cap	30x10 tubes	010-501-564
	100mL <i>DigiTUBE</i> s with RackLock, blue cap	30x10 tubes	010-501-561
	100mL <i>DigiTUBE</i> s with RackLock, clear cap	30x10 tubes	010-501-562



Disposable Watch Glasses

Description	Quantity	Part Number
Disposable Polypropylene Watch Glasses, 50 mL	1000	010-500-081
Disposable Polypropylene Watch Glasses, 100 mL	500	010-501-081

Promote refluxing or protect your samples from environmental contamination during digestions with our high-quality, disposable watch glasses. These products are available to pair with 50 mL or 100 mL disposable *DigiTUBE*s. Our watch glasses are manufactured using virgin polypropylene and designed for easy handling.

Each watch glass is equipped with a ribbed step to keep it securely in place during use and a handle to make placement and removal easier while also preventing cross-contamination during use.



Storage Racks

Description	Quantity	Part Number
Plastic Storage Racks for 50 mL <i>DigiTUBE</i> s	5	010-500-026
Foam Storage Racks for 100 mL <i>DigiTUBE</i> s	5	010-501-028

