

PremiSolv

ICP Solvent

PremiSolv ICP solvent is ideal for diluting oil and organic liquid samples for analysis by ICP-OES, RDE-OES or UVF.

Analyzed by ICP for any trace presence of sulfur and 40 metals, PremiSolv ICP solvent delivers the reliable performance you need. In addition, PremiSolv ICP solvent is extremely low in toxicity and odor compared to alternative industrial solvents such as kerosene, increasing lab safety.

PremiSolv ICP solvent provides unparalleled overall solution stability and convenience as a result of:

- ▶ Proven oil-based fluid and hydrocarbon solubility
- ▶ Robust micro-filtration of every lot, ensuring cleanliness
- ▶ Convenient pack sizes and easy-to-use containers
- ▶ Simple shipping due to non-hazardous classification
- ▶ Pre-mixed PremiSolv solvents with internal standard (IS) cobalt at 40 ppm are also available

Whether you are performing standard methods such as ASTM D5185 or D6595, or preparing routine calibration standards, PremiSolv meets your total requirements.



Description	Size	Part number
PremiSolv	1 Drum (208L)	150-700-001
PremiSolv	5 Gallons (18.9 L)	150-700-002
PremiSolv	1 Gallon (3.78 L)	150-700-003
PremiSolv with IS (Co 40 ppm)	1 Drum (208L)	150-005-274
PremiSolv with IS (Co 40 ppm)	5 Gallons (18.9 L)	150-005-272
PremiSolv with IS (Co 40 ppm)	1 Gallon (3.78 L)	150-005-271

