

Plate Specifications

for AnalytiChem UK

Sizes and Dimensions

The **Redipor**® range of plated media from **AnalytiChem UK** (formerly Cherwell Laboratories) is available in 140 mm, 90 mm and 55 mm plates. All dishes are manufactured from clear polystyrene and are vented.



Our range of Petri dishes are provided with a standard fill volume. However, these are customizable to have deeper fill volumes. **AnalytiChem UK** also has the ability to provide bespoke packaging solutions. Please contact us for further information, to discuss your requirements.

Technical Data

Description

Type: Petri dish

Material / Type: Clear polystyrene

Nominal Size	No. of Vents	Maximum Fill	Overall Dimensions (empty, incl. ports)
55 mm	Six (6)	9.5 mL	14.2 (H) x 58.0 (Dia) mm
90 mm	Single (standard) or Triple (optional)	32 mL	14.2 (H) x 92.0 (Dia) mm
140 mm	Triple	130 mL	21.4 (H) x 143.5 (Dia) mm

Type: Contact plates **Material / Type:** Clear polystyrene with grip lid

Nominal Size	No. of Vents	Maximum Fill	Overall Dimensions (empty, incl. ports)
55 mm	Triple	17 mL	14.2 (H) x 65.8 (Dia) mm

Redipor® TwistLock plates are also available. Our range of plated prepared media with locking lids, designed to enhance sample security, data integrity and workflow efficiency. Contact us for further details.

Terms & Conditions

All specifications are approximate, where precise data is required please enquire about samples.

We reserve the right to amend specifications without any prior notice. E&OE.

