redipor
an analytichem brand

Universal Operator Broth Transfer Kit



A Kit to Assess Operator Technique

The **Redipor®** Operator Broth Test kit is a prepackaged set of components created in accordance with guidance from the UK Pharmaceutical Aseptic Services Committee.

The guidance was put together by *Tim Sizer, Robert Duncombe, Richard Needle, Graham Sewell*; the purpose of the kit is to aid in the assessment of operators undertaking aseptic preparations, to maintain the sterility of materials during preparation and to ensure that minimal bio-burden is introduced.

Applications

The Redipor Universal kit is for use in high level clean room environments as a standard assessment of an operator's ability to maintain the sterility of materials during the preparation of aseptically prepared injectable dose forms, including to form part of routine operator monitoring.

redipor
an analytichem brand

Features & Benefits

Each kit contains a 50 mL Tryptic Soy Broth injection vial, three empty 30 mL sterile injection vials, a 100 mL Tryptone Soya Broth bag and a 10 mL Tryptone Soya Broth ampule.

The EVA infusion bag is particularly appropriate to TPN compounding. All components are supplied by AnalytiChem and packaged into a purposefully designed corrugated plastic box for ease of use.

As with all our Redipor prepared media, all batches are accompanied with detailed Certificate of Analysis.

Technical Data and Ordering Information

Qty.	Kit contents
1	Tryptic Soy Broth 50 mL Injection vial
3	Empty 30 mL Injection vials (sterile)

Qty.	Kit contents
1	Tryptone Soya Broth 100 mL broth bag
1	10 mL Tryptone Soya Broth Ampoule

Part Number	Description	Qty.
4.0286	Redipor Universal Operator Broth Transfer Kit	Each

Product Specifications

Packaging	All items packed in a corrugated plastic case
Dimensions (mm)	176 x 132 x 137 (L x W x H)
Weight (complete kit; kg)	0.5

Terms & Conditions

We reserve the right to amend specifications without notice. E&OE

AnalytiChem UK Ltd.
sales.uk@analytichem.com
+44 (0) 1869 355500

7 & 8 Launton Business Centre,
Murdock Road, Bicester,
OX26 4XB, UK

© 2026 AnalytiChem. All rights reserved.
All trademarks and registered trademarks are the property of their respective owners.
www.analytichem.com

AC_LIT_UOBTK_UK-A4-EN-001